

# LabVIEW Report

## 28/01/2021

P. Hacker

**HELMHOLTZ**  
SPITZENFORSCHUNG FÜR  
GROSSE HERAUSFORDERUNGEN



This work has been carried out within the framework of the EUROfusion Consortium and has received funding from the Euratom research and training programme 2014-2018 and 2019-2020 under grant agreement No 633053. The views and opinions expressed herein do not necessarily reflect those of the European Commission.

- reviewed the LabVIEW DAQ routine to minimize lag, calculation overhead and optimize comfort of use for simplicity
- started from 'original' master branch that included all changes from frequency selection, power calculation, register settings and real time feedback
- build two new DAQ programs that include a number of changes:
  - removed all disable and if structures that were originally used for debugging
  - only manual acquisition time input
  - two trigger settings: start of VI time, CODAC last trigger
  - timing inputs to [s], time display as needed in [s, ms, ns]
  - samples to exclude from acquisition fix to 1000, included in start/time calculation of acquisition itself
  - DAC range settings now only either global or individual via file input (10-80mV)

- reviewed the LabVIEW DAQ routine to minimize lag, calculation overhead and optimize comfort of use for simplicity
- started from 'original' master branch that included all changes from frequency selection, power calculation, register settings and real time feedback
- build two new DAQ programs that include a number of changes:
  - calibration values for integrity check can be calculated from current measurement or either input for lab tests via file (provided one with OP1.2b 95% quantile)
  - calibration values now in [ $A^2$ , Ohm, s]
  - raw signal display (only, no adjusted or power per channel) in [mV]
  - offset measurement for integrity check, display in [mV]
  - reduced archive upload, only necessary data

# Reviewed LabVIEW DAQ routine

- reviewed the LabVIEW DAQ routine to minimize lag, calculation overhead and optimize comfort of use for simplicity
- started from 'original' master branch that included all changes from frequency selection, power calculation, register settings and real time feedback
- build two new DAQ programs that include a number of changes:

- new timing settings: assuming that there is an intrinsic timing error from internal measurement, i.e. inside DAQ loop, of  $\sim 0.04\text{ms}$

HEX	DECIMAL	FREQ [Hz]	T [ms]
=====			
x4F	79	1318	0.7584
x63	99	1052	0.9504
xA2	162	643.0	1.5552
x148	328	317.6	3.1488
x294	660	157.8	6.336

## Raw wRTF

- real time feedback plot with full  $P_{\text{rad}}$  of both camera arrays and single channel/fast  $P_{\text{rad}}$
- geometry input via file
- power calculation now only via FIFO, no Savitzky Golay

## Raw w/oRTF

- removed voltage out, power calculation,  $P_{\text{rad}}$  routine, geometry input entirely




# Reviewed LabVIEW DAQ routine: Performance

master

Raw\_wRTF


Raw\_w/oRTF

0.8ms




# Reviewed LabVIEW DAQ routine: Performance

master




Raw\_wRTF

1.0ms




Raw\_w/oRTF




# Reviewed LabVIEW DAQ routine: Performance

master




Raw\_wRTF

1.6ms




Raw\_w/oRTF






# Reviewed LabVIEW DAQ routine: Performance

master




Raw\_wRTF

3.2ms




Raw\_w/oRTF




# Reviewed LabVIEW DAQ routine: Performance

master




Raw\_wRTF

6.4ms




Raw\_w/oRTF



# Reviewed LabVIEW DAQ routine: Performance


unloaded



raw\_woRTF

0.8ms

loaded




# Reviewed LabVIEW DAQ routine: Performance


raw\_woRTF

1.0ms

unloaded




loaded



# Reviewed LabVIEW DAQ routine: Performance


samples per laser half period

master




Raw\_wRTF

0.8ms




Raw\_w/oRTF



# Reviewed LabVIEW DAQ routine: Performance


samples per laser half period

master




Raw\_wRTF

1.0ms




Raw\_w/oRTF



# Reviewed LabVIEW DAQ routine: Performance


## samples per laser half period

master




Raw\_wRTF

1.6ms



Raw\_w/oRTF






# Reviewed LabVIEW DAQ routine: Performance


samples per laser half period

master




Raw\_wRTF

3.2ms



Raw\_w/oRTF





# Reviewed LabVIEW DAQ routine: Performance




samples per laser half period

master

Raw\_wRTF

Raw\_w/oRTF

6.4ms




# Reviewed LabVIEW DAQ routine: Performance


unloaded

raw\_woRTF

0.8ms



loaded





# Reviewed LabVIEW DAQ routine: Performance

unloaded

raw\_woRTF

1.0ms

loaded



- [https://git.ipp-hgw.mpg.de/pih/LabVIEW\\_QSB](https://git.ipp-hgw.mpg.de/pih/LabVIEW_QSB) holds all relevant changes and versions under branches, i.e. master, raw\_wRTF and raw\_w/oRTF
- Changelog.md includes all information about the VI I was able to acquire and learn throughout the optimization process
- versions as of today (02/05/2020) also located in QSB\_Bolometry E5 backup server