



Max-Planck-Institut für Plasmaphysik

Report 04/16/2020

P. Hacker









This work has been carried out within the framework of the EUROfusion Consortium and has received funding from the Euratom research and training programme 2014-2018 and 2019-2020 under grant agreement No 633053. The views and opinions expressed herein do not necessarily reflect those of the European Commission.

Setup







rectangular aperture at location of old one (HBCm) with plane normal exactly pointing at (0.0, 0.0, 0.0)

detector array equidistant to aperture and to each other behind pinhole and inline with plane normal

detector plane normal has no angle with pointing vector between aperture and center









Setup





➤ Lines of sight from this configuration are symmetric (top-bottom) down to the numerical accuracy of a 53 bit encoded double precision float






sum of opposing channels (0-31, 1-30,...) in (x,y,z)-space

Emissivity on Fluxsurface Cells








- right almost entirely symmetric emissivity distribution on fluxsurfaces, except artifact in front of aperture
- > slightly wider opening angle on line of sight-fan

Forward Integral: Standard Case vs. Artificial Array







- ris located not exactly at 108 degrees, but rather 107.9)
- > geometrically very different to 'standard' case
- >more peaked towards edge
- > far inside the LCFS even for low radiation fraction cases

2D Tomography: First Results






geometry files are now being read correctly

➤ anisotropy factors set to 1.5 for divertor (?) or island and 0.8 for core-like regions

 \geq t = 3.421s <-> frad=100%




2D Tomography: First Results





>anisotropy factors set to 1.0 and 1.0

> again: t = 3.421s (frad=100%)



Agenda





Not in order:

>start to finish RSI paper (for internal reviewing)

>understand anisotropy factor, hence structure of tomogram

>create and feed phantom radiation profiles to inversion method using that knowledge



1.: 'Fix' Bolometer Planar Error






> guided by lower half of detector array
(green)


➤ observe difference (red) between opposite channel (e.g. CH#0 <-> CH#31) through rotating around aperture normal by 180° and measuring angle

> transforming second channel through rotating it by angle from before (blue)

➤ only really easy for HBCm, because central aperture axis alignment; VBC cameras not possible







2.: Toroidal Transformation to Axis Plane






➤ guide is normal of plane constructed by detector fan (blue)

transforming all channels so that normal points in/at toroidal direction (red)

➤ transformation for each channel individual (see previous argument (1.)) (orange)

>done for all cameras individually






3.: Tilting the Detector Fan Up/Down





≻ take results of (2.) and tilt the entire fan including the aperture poloidally (grey to blue and black to red)




change included)

Centering of Aperture and Random Error







>center the aperture center (black) on z = 0 axis and shift entire array accordingly to red



> degree of tilt 0.5° - 2.0°, but omnidirectional





-1.3193983963943924 834 834 834 838 839 849 841 841 84